



SINUSCOPY
-CATALOG-



Our new handles with axial cleaning port and suction port

Model II with cleaning port



Model Suction II with suction port



Our new handle with axial cleaning port and 360° rotatable





INNOVATION handle design
45° – cleaning port

Besides our existing rinsable instruments with axial LUER-LOCK connection, from now one we can offer also another variation of a well-designed 45° – LUER-LOCK connector,
With also excellent hygienic cleaning effect.



Blakesley forceps through-cutting, jaw straight



INSERT DETAILS		Art. No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	2,5 mm (Size 0)	SIN10-0A	SIN10-0B	SIN10-0G
130 mm	3,0 mm (Size 1)	SIN10-1A	SIN10-1B	SIN10-1G
130 mm	3,6 mm (Size 2)	SIN10-2A	SIN10-2B	SIN10-2G
130 mm	4,2 mm (Size 3)	SIN10-3A	SIN10-3B	SIN10-3G

All models with integrated cleaning port (Luer lock connection)



Blakesley forceps through-cutting, jaw 45° upwards



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	2,5 mm (Size 0)	SIN20-0A	SIN20-0B	SIN20-0G
130 mm	3,0 mm (Size 1)	SIN20-1A	SIN20-1B	SIN20-1G
130 mm	3,6 mm (Size 2)	SIN20-2A	SIN20-2B	SIN20-2G
130 mm	4,2 mm (Size 3)	SIN20-3A	SIN20-3B	SIN20-3G

All models with integrated cleaning port (Luer-Lock connection)



Grünwald-Henke forceps through-cutting, jaw straight



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	3,0 mm (Size 1)	SIN30-1A	SIN30-1B	SIN30-1G
130 mm	3,6 mm (Size 2)	SIN30-2A	SIN30-2B	SIN30-2G

All models with integrated cleaning port (Luer Lock connection)



SUCTION Blakesley through- cutting forceps



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model Suction I 	Model Suction II 	
130 mm	3,0 mm (Size 1)	SIN50-1C	SIN50-1D	
130 mm	3,6 mm (Size 2)	SIN50-2C	SIN50-2D	

All models with integrated suction channel(Luer-Lock connection)



Binner polyp forceps



INSERT DETAILS		Art.No.		
Length (in mm)	Shape of shaft	Model I 	Model Rotatable II 	
130 mm	Straight	SIN80-0A	SIN80-0H	
130 mm	45° up curved	SIN80-45A		
130 mm	80° up curved	SIN80-80A		

All models with integrated cleaning port (Luer-Lock connection)



Biopsy & Grasping forceps with Y-handle, malleable shaft



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model III 	Model VI 	
145 mm	3,0 mm	SIN90-3E	SIN90-3F	
145 mm	4,0 mm	SIN90-4E	SIN90-4F	

All models with integrated cleaning port (Luer-Lock connection)



SUCTION - Blakesley forceps straight with integrated suction channel



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model Suction I 	Model Suction II 	
130 mm	3,0 mm (Size 1)	SIN60-1C	SIN60-1D	
130 mm	3,6 mm (Size 2)	SIN60-2C	SIN60-2D	
130 mm	4,2 mm (Size 3)	SIN60-3C	SIN60-3D	

All models with integrated suction channe(Luer-Lock connection)



Blakesley forceps straight



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	2,5 mm (Size 0)	SIN100-0A	SIN100-0B	SIN100-0G
130 mm	3,0 mm (Size 1)	SIN100-1A	SIN100-1B	SIN100-1G
130 mm	3,6 mm (Size 2)	SIN100-2A	SIN100-2B	SIN100-2G
130 mm	4,2 mm (Size 3)	SIN100-3A	SIN100-3B	SIN100-3G
130 mm	4,8 mm (Size 4)	SIN100-4A	SIN100-4B	SIN100-4G

All models with integrated cleaning port (Luer-Lock connection)



Blakesley-Wilde forceps 45° upwards , 90° right horizontal to position of handle



INSERT DETAILS		Art.No		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
130 mm	3,6 mm (Size 1)	SIN112-1A	SIN112-1B	

All models with integrated cleaning port (Luer-Lock connection)



Blakesley forceps jaw straight, shaft 25° upwards curved



INSERT DETAILS		Art.No		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
130 mm	2,5 mm (Size 0)	SIN112-1A	SIN112-1B	
130 mm	3,0 mm (Size 1)	SIN105-1A	SIN105-1B	
130 mm	3,6 mm (Size 2)	SIN105-2A	SIN105-2B	

All models with integrated cleaning port (Luer-Lock connection)



Blakesley-Wilde forceps 45° upwards



INSERT DETAILS		Art.No		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	2,5 mm (Size 0)	SIN110-0A	SIN110-0B	SIN110-0G
130 mm	3,0 mm (Size 1)	SIN110-1A	SIN110-1B	SIN110-1G
130 mm	3,6 mm (Size 2)	SIN110-2A	SIN110-2B	SIN110-2G
130 mm	4,2 mm (Size 3)	SIN110-3A	SIN110-3B	SIN110-3G

All models with integrated cleaning port (Luer-Lock connection)



Blakesley-Wilde forceps 45° upwards, shaft 25° curved up



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
130 mm	2,5 mm (Size 0)	SIN115-0A	SIN115-0B	
130 mm	3,0 mm (Size 1)	SIN115-1A	SIN115-1B	
130 mm	3,6 mm (Size 2)	SIN115-2A	SIN115-2B	

All models with integrated cleaning port (Luer-Lock connection)



Blakesley-Wilde forceps, jaw 90° upwards



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	3,6 mm (Size 1)	SIN120-1A	SIN120-1B	SIN120-1G
130 mm	4,2 mm (Size 2)	SIN120-2A	SIN120-2B	SIN120-2G

All models with integrated cleaning port (Luer-Lock connection)



Circular cutting punch egg-form with Y-Handle



STRAIGHT SHAFT

INSERT DETAILS		Art.No.
Length (in mm)	Diameter of jaw (in mm)	Model III 
180 mm	3,5 mm	SIN130-0E
180 mm	4,5 mm	SIN135-0E

All models with integrated cleaning port (Luer-Lock connection)

65° CURVED UP IN SHAFT

INSERT DETAILS		Art.No.
Length (in mm)	Diameter of jaw (in mm)	Model III 
180 mm	3,5 mm	SIN130-65E
180 mm	4,5 mm	SIN135-65E

All models with integrated cleaning port (Luer-Lock connection)



Circular cutting punch with window and Y-Handle



STRAIGHT SHAFT

INSERT DETAILS		Art.No.
Length (in mm)	Diameter of jaw (in mm)	Model III 
180 mm	3,5 mm	SIN390-0E
180 mm	4,5 mm	SIN395-0E

All models with integrated cleaning port (Luer-Lock connection)

65° CURVED UP IN SHAFT

INSERT DETAILS		Art.No.
Length (in mm)	Diameter of jaw (in mm)	Model III 
180 mm	3,5 mm	SIN390-65E
180 mm	4,5 mm	SIN395-65E

All models with integrated cleaning port (Luer-Lock connection)



Takahashi forceps



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
180 mm	2,5 x 8 mm	S140-0A	S140-0B	S140-0G
180 mm	2,5 x 8 mm	S141-0A	S141-0B	S141-0G
180 mm	2,5 x 8 mm jaw 30°	S142-30A	S142-30B	S142-30G

All models with integrated cleaning port (Luer-Lock connection)



SUCTION - Takahashi forceps with integrated suction channel



INSERT DETAILS		Art.No		
Length (in mm)	Diameter of jaw (in mm)	Model Suction I 	Model Suction II 	
130 mm	2,5 x 8 mm	S70-0C	S70-0D	
180 mm	2,5 x 8 mm	S71-0C	S71-0D	

All models with integrated suction channel(Luer-Lock connection)



Heuwieser antrum polyp grasping forceps



INSERT DETAILS		Art.No.	
Length (in mm)	Shape	Model VI 	Model Rotatable II 
100 mm	90° curved 120° backward opening	SIN160-0F	SIN160-0H
100 mm	120° curved 150° backward opening	SIN161-0F	SIN161-0H

All models with integrated cleaning port (Luer-lock connection)



Hilyard post biopsy forceps, round jaw



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
140 mm	2,5 mm	SIN170-25A	SIN170-25B	SIN170-25G
130 mm	3,0 mm	SIN171-30A	SIN171-30B	SIN171-30G

All models with integrated cleaning port (Luer-Lock connection)



Watson-William etmoid forceps



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
130 mm	3,0 mm	SIN180-0A	SIN180-0B	

All models with integrated cleaning port (Luer-Lock connection)



Hartmann through-cutting forceps



INSERT DETAILS		Art.No		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
130 mm	5,0 mm	SIN200-5A	SIN200-5B	
130 mm	7,0 mm	SIN200-7A	SIN200-7B	

All models with integrated cleaning port (Luer-Lock connection)



Schmeden through-cutting forceps, triangular jaw



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
130 mm	Size 1	SIN210-1A	SIN210-1B	
130 mm	Size 2	SIN210-2A	SIN210-2B	

All models with integrated cleaning port (Luer-Lock connection)



Kerrison punch 90° up



INSERT DETAILS

Art.No.

INSERT DETAILS		Art.No.
Length (in mm)	Diameter of jaw (in mm)	
180 mm	2,0 mm	SIN220-1
180 mm	3,0 mm	SIN220-2
180 mm	4,0 mm	SIN220-3



Nasal scissors



INSERT DETAILS		Art.No.		
Length (in mm)	Shape	Model IV 	Model I 	
130 mm	Straight, serrated	SIN230-SF	SIN230-SA	
130 mm	Left curved, serrated	SIN230-LF	SIN230-LA	
130 mm	Right curved, serrated	SIN230-RF	SIN230-RA	

All models with integrated cleaning port (Luer-Lock connection)



Frontal and Sphenoidal, through-cutting, double action punch



INSERT DETAILS		Art.No.		
Length (in mm)	Bite of jaw (in mm)	Model V 	Model IV 	
140 mm	2,0 mm vertical	SIN240-55I	SIN240-55F	
140 mm	2,0 mm vertical	SIN240-90I	SIN240-90F	
140 mm	2,0 mm horizontal	SIN241-55I	SIN241-55F	
140 mm	2,0 mm horizontal	SIN241-90I	SIN241-90F	

All models with integrated cleaning port (Luer-Lock connection)



Maxillary forceps, double action jaw



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model V 	Model IV 	
150 mm	2,5 mm straight	SIN250-SI	SIN250-SF	
120 mm	3,0 mm 70° down vertical	SIN251-70I	SIN251-70F	
110 mm	3,0 mm 110° down horizontal	SIN252-110I	SIN252-110F	

All models with integrated cleaning port (Luer-Lock connection)



Giraffe - Frontal sinus forceps, double action jaws



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model V 	Model IV 	
140 mm	4,0 mm 55° curved up, vertical	SIN260-4-55I	SIN260-4-55F	
140 mm	4,0 mm 55° curved up, horizontal	SIN261-4-55I	SIN261-4-55F	
140 mm	3,0 mm 70° curved up, vertical	SIN260-3-70I	SIN260-3-70F	
140 mm	3,0 mm 70° curved up, horizontal	SIN261-3-70I	SIN261-3-70F	
140 mm	3,0 mm 110° curved up, vertical	SIN260-3-110I	SIN260-3-110F	
140 mm	3,0 mm 110° curved up, horizontal	SIN261-3-110I	SIN261-3-110F	
140 mm	2,0 mm 90° curved up, vertical	SIN260-2-90I	SIN260-2-90F	
140 mm	2,0 mm 90° curved up, horizontal	SIN261-2-90I	SIN261-2-90F	

All models with integrated cleaning port (Luer-Lock connection)



Retro-Antrum punch, 360 ° rotatable



INSERT DETAILS		Art.No.		
Length (in mm)	Bite of jaw (in mm)	Model Rotatable III	Model Rotatable I	
100 mm	1,5 mm			
100 mm	2,5 mm	SIN280-15J SIN281-25J	SIN280-15G SIN281-25G	

All models with integrated cleaning port (Luer-Lock connection)



TOBEY Retro-Antrum punch, 360° rotatable, shaft slightly curved



INSERT DETAILS		Art.No.		
Length (in mm)	Bite of jaw (in mm)	Model Rotatable III	Model Rotatable I	
				
100 mm	2,0 mm	SIN288J	SIN288G	

All models with integrated cleaning port (Luer-Lock connection)



Antrum-Backbiter punch



INSERT DETAILS		Art.No.		
Length (in mm)	Diameter of jaw (in mm)	Model I 	Model II 	
100 mm	1,5 mm straight	SIN290-15A	SIN290-15B	
100 mm	1,5 mm left	SIN291-15A	SIN291-15B	
100 mm	1,5 mm right	SIN292-15A	SIN292-15B	
100 mm	2,5 mm straight	SIN293-25A	SIN293-25B	
100 mm	2,5 mm left	SIN294-25A	SIN294-25B	
100 mm	2,5 mm right	SIN295-25A	SIN295-25B	

All models with integrated cleaning port (Luer-Lock connection)



Struycken forceps



INSERT DETAILS		Art.No		
Length (in mm)	Bite of jaw (in mm)	Model I 	Model II 	Model Rotatable I 
130 mm	1,5 x 9 mm	SIN300-1A	SIN300-1B	SIN300-1G
180 mm	1,5 x 9 mm	SIN301-1A	SIN301-1B	SIN301-1G

All models with integrated cleaning port (Luer-Lock connection)



SCHWAGER through-cutting punch for Orbitaboden



INSERT DETAILS

Art.No.

INSERT DETAILS		Art.No.	
Length (in mm)	Bite of jaw (in mm)	Model I 	Model II 
130 mm	2,0 mm 90° left in jaw	SIN310-LA	SIN310-LB
130 mm	2,0 mm 90° right in jaw	SIN311-RA	SIN311-RB
140 mm	2,0 mm 90° shaft curved up	SIN312-UA	SIN312-UB

All models with integrated cleaning port (Luer-Lock connection)



SUCTION - Blakesley forceps, jaw straight & 45°up



with separately suction channel

INSERT DETAILS		Art.No.
Length (in mm)	Diameter of jaw (in mm)	
130 mm	3,6 mm (Size 1) straight	SIN330-1
130 mm	4,2 mm (Size 2) straight	SIN330-2
130 mm	3,6 mm (Size 1) 45° up	SIN340-1
130 mm	4,2 mm (Size 2) 45° up	SIN340-2

All models with integrated cleaning port (Luer-Lock connection)



Flexible instruments for Nasopharyngoscope



Alligator forceps



Biopsy forceps

INSERT DETAILS

Art.No.

Length (in mm)	Diameter of jaw (in mm)	Model VI	Model VII	
				
400 mm	1,0 mm Biopsy forceps	SIN360-BK	SIN360-BL	
300 mm	1,2 mm Biopsy forceps	SIN362-BK	SIN362-BL	
400 mm	1,0 mm Alligator forceps	SIN363-AK	SIN363-AL	
300 mm	1,2 mm Alligator forceps	SIN361-AK	SIN361-AL	

Model VI with integrated cleaning port (Luer-Lock connection)



Flexible instruments for Nasopharyngoscope, demountable



Alligator forceps



Biopsy forceps

INSERT DETAILS

Art.No.

Length (in mm)	Diameter of jaw (in mm)	Model VIII	Insert only		Handle only	
						
400 mm	1,0 mm Biopsy forceps	SIN365-BM	SIN365-B		SIN-CWG	
300 mm	1,2 mm Biopsy forceps	SIN367-BM	SIN367-B		SIN-CWG	
400 mm	1,0 mm Alligator forceps	SIN368-AM	SIN368-A		SIN-CWG	
300 mm	1,2 mm Alligator forceps	SIN366-AM	SIN366-A		SIN-CWG	

System 1



System 2



INNOVATION

- Large variety of interchangeable forceps and scissors inserts.
- All types of instrument inserts and handles are compatible. - Assembly and disassembly at the push of a button. Safe fixation of the instrument inserts by click mechanism.
- ergonomic handle
- Easy and efficient cleaning/sterilisation by axial irrigation.
- Re-useable, evident cost reduction compared to standard non-detachable instruments.
- Self-locking spring, preventing unintentional opening of the jaws.
- Irrigation of handle mechanism by luer-lock cone.
- innovative patented system

INFORMATIONS ABOUT THE COMPANY:

ENDO-PASSION is your competent partner for innovative products in the field of flexible endoscopy. The force of innovation is drawn from years of experience in production of endoscopic instruments.

In our headquarters in Germany, we produce instruments using advanced production processes. We use only materials, manufacturing and testing processes, which can withstand the highest quality requirements.

Our experienced staff develop and manufacture with the highest precision innovative and user-friendly solutions for endoscopy



ENDO-PASSION GmbH
Faulenbachstr. 15, 78573 Wurmlingen
Tel. 07461-1617771, Fax 07461-1617772
E-Mail: info@endo-passion.de
www.endo-passion.de